

TRANSPORT MEDIA

COMPATIBILITY



Table evaluates compatibility of transport media with Fluxergy's COVID-19 Test Kit.

Transport Media Name	Catalog Number	Media Compatible with Reaction Mixture and Test Card	Sold in 3mL	Glass Beads Present	Includes Stability Agents for Viral Storage
Copan UTM 3mL™		✓	✓	✓	✓
Puritan UT-300™	UT-300	✓	✓	✓	✓
BD UVT 3 mL™	220527, 220528, 220529, 220531	✓	✓	✓	✓
Remel MicroTest™ M4RT™	R12552	✓	✓	✓	✓
Citoswab VTM™	780028	✓	✓		✓
Medschenker STM™	STM 30-A	✓	✓		✓
Avantik VTM™	GL4692	✓	✓		✓
Hardy	3C036NHL, 3C037NHL, 3C038NHL	TBD	✓		✓
CDC	SOP# DSR-052-05	TBD	✓		✓
PBS without Magnesium or Calcium		TBD	✓		

Fluxergy's products comply with DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices. The devices are not endorsed or approved outside of the European Economic Area.

Version Date: February 19, 2021

Fluxergy
30 Fairbanks, Suite 110
Irvine, CA 92618 - USA
(949) 305-4201
customersupport@fluxergy.com
www.fluxergy.com