

FLUXERGY MATERIALS CHECKLIST

COVID-19 TEST KIT



Provided Materials:	Included
Fluxergy Analyzer	Contains Fluxergy Analyzer, 1 Power Cord and 1 ethernet cord
Fluxergy COVID-19 Reaction Mixture	Contains 100 frozen vials of single-use COVID-19 Reaction Mixture
Fluxergy COVID-19 Test Card	Contains 100 single-use COVID-19 Test Cards
Fluxergy Works Software	Registration key to create access logins for multiple users

*Material specifics can be found in [Instructions for Use](#).

Materials **NOT** Provided*:

Fluxergy Analyzer Platform:
PC with Fluxergy Works Installed
Network or Router to connect Analyzer and PC
Ethernet Cables (Quantity for 1 Analyzer and PC)
2D Laser Barcode Scanner
Power Strip (Quantity for 1 Analyzer, PC and Workflow Materials)
Testing Material:
VTM (3 mL) and swabs for Sample Collection. <i>Please see our List of Accepted VTMs.</i>
Quality Controls. <i>Please see our List of Accepted Quality Controls.</i>
Laboratory Supplies:
200uL calibrated pipette and associated filtered tips
20 or 30uL calibrated pipette and associated filtered tips
Mini Vortex Mixer
Mini Centrifuge (2000g)
Manual defrost Freezer (-20°C) to store reaction mixtures
Refrigerator for potential storage of samples
PPE and cleaning solutions for safe testing
Biosafety Cabinet (if desired)

*Material specifics can be found in [Instructions for Use](#).



*Only Fluxergy Test Kit COVID-19 has been approved for CE-IVD and is available for any markets that accept CE marking as the valid regulatory approval. All other products or any other markets are For Research Use Only (RUO) and are not for use in diagnostic procedures.