


FLUXERGY MATERIALS CHECKLIST

EQUINE HERPES VIRUS – 1 (EHV-1)



Provided Materials:	Included
Fluxergy Analyzer	Contains Fluxergy Analyzer, 1 Power Cord and 1 ethernet cord
Fluxergy EHV-1 Reaction Mix	Contains 10 frozen vials of EHV-1 Reaction Mix
Viral Transport Media (3 mL)	Contains 10 vials of Viral Transport Media (3 mL)
Sterile Flock Swabs	Contains 10 individually wrapped Flock Swabs
Fluxergy EHV-1 Test Card	Contains 10 EHV-1 Test Cards
Fluxergy Works Software	Registration key to create access logins for multiple users

***Material specifics can be found in Instructions for Use.*

Materials NOT Provided*:

Fluxergy Analyzer Platform:
PC with Fluxergy Works Installed**
Network or Router to connect Analyzer and PC **
Ethernet Cables (Quantity for 1 Analyzer and PC) **
2D Laser Barcode Scanner **
Power Strip (Quantity for 1 Analyzer, PC and Workflow Materials)
Testing Material:
Quality Controls. ** <i>Material specifics can be found in Instructions for Use.</i>
Laboratory Supplies:
200uL calibrated pipette and associated filtered tips **
Mini Vortex Mixer **
Mini Centrifuge (2000g) **
Manual defrost Freezer (-20°C) to store reaction mixtures
Refrigerator for potential storage of samples
PPE and cleaning solutions for safe testing
Biosafety Cabinet (if desired)

***Material specifics can be found in Instructions for Use.*



No Fluxergy products have been cleared or approved by the USDA for veterinary diagnostic use. Fluxergy's current veterinary product offerings are limited to individual laboratory components for use by experienced professional users in developing their own laboratory tests for research or animal use.